



## Buffer Powder Pillows, pH 6.86, pk/15

**Product #:** 1409895  
**ZAR Price:** Contact Hach  
**Available**

### pH Buffer standard solutions for laboratory and field analyses

pH Buffer Powder Pillows for pH electrode calibration. Each powder Pillow prepares 50 mL of fresh buffer solution. Colourless. pH is 6.86  $\pm 0.02$  at 25°C. Pack of 15 Powder Pillows. Pack of 15 Powder Pillows

#### Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

### Specifications

|                     |                                     |
|---------------------|-------------------------------------|
| Color coded?:       | No                                  |
| Description:        | pH Buffer Calibration reagent       |
| NIST Traceable:     | Yes                                 |
| Package Type:       | Powder Pillows                      |
| pH Value(s):        | 6.86 $\pm 0.02$ pH at 25°C          |
| Quantity:           | 15 pcs                              |
| Special Remarks:    | Discard prepared solution after use |
| Storage conditions: | 10 - 25 °C (protect from light)     |