



## Buffer Powder Pillows, pH 8.30, 50 mL, pk/25

**Product #:** 89868  
**ZAR Price:** Contact Hach  
**Ships within 3-5 weeks**

Buffer for preparing monochloramine standards and detecting visual alkalinity endpoints. pH 8.30, Sample Size: 50 mL sample. Pack of 25 Powder Pillows.

### Fresh and reliable

Unit-dose powder pillows assure buffer freshness and eliminate the possibility of contamination.

---

### Specifications

Concentration:	pH 8.3
Description:	Buffer Powder Pillows, 25/pk
Platform:	Powder Pillows