



## Buffer Solution, pH 4.01, Colour-coded Red, 500 mL

**Product #:** 2283449  
**ZAR Price:** Contact Hach  
**Available**

### pH Buffer standard solutions for laboratory and field analyses

pH Buffer Standard Solution for pH electrode calibration. Colour-coded pH Buffer solution for easy identification. pH is  $4.01 \pm 0.02$  at  $25^{\circ}\text{C}$ . 500 mL bottle

#### Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

### Specifications

|                     |  |
|---------------------|--|
| Color coded?:       | Yes  |
| Description:        | pH Buffer Calibration reagent  |
| NIST Traceable:     | Yes  |
| Package Type:       | Bottle   |
| pH Value(s):        | $4.01 \pm 0.02$ pH at $25^{\circ}\text{C}$   |
| Quantity:           | 500 mL   |
| Special Remarks:    | Observed shelf-life after opening, following recommended storage and use: 6 months |
| Storage conditions: | $10 - 25^{\circ}\text{C}$ (protect from light)                                     |