



## Singlet Single-Use pH Buffer Kit, pH 4.01 & 7.00, pk/2x10

**Product #:** 2769920  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

### Singlet Solution Packs

pH Buffer Standard Solution for pH electrode calibration. Sealed foil pouches of pH buffers. pH is  $\pm 0.02$  at 25°C. Pack of 10 of pH 4.01 and 10 of pH 7.00, 20 mL each

#### Ideal for on-site analyses

Sealed pouches are easy to transport - no risk of spillage.

#### Easy to Use

Just open the pouch and insert the electrode.

#### Fresh and reliable

Foil pouches prevent contamination and air exposure.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

---

### Specifications

Color coded?:	Yes
Description:	pH Buffer Calibration reagent
NIST Traceable:	Yes
Package Type:	Singlet
pH Value(s):	4.01 & 7.00 $\pm 0.02$ pH at 25°C
Quantity:	2 pks of 10 pcs
Special Remarks:	Discard solution after use
Storage conditions:	10 - 25 °C (protect from light)