



Singlet Single-Use pH Buffer, pH 7.00, pk/20

Product #: 2770120
ZAR Price: Contact Hach
Ships within 3 days

Singlet Solution Packs

pH Buffer Standard Solution for pH electrode calibration. Sealed foil pouches of pH buffers. pH is ± 0.02 at 25°C. Pack of 20 of pH 7.00, 20 mL each

Ideal for on-site analyses

Sealed pouches are easy to transport - no risk of spillage.

Easy to Use

Just open the pouch and insert the electrode.

Fresh and reliable

Foil pouches prevent contamination and air exposure.

Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

Specifications

Color coded?:	Yes
Description:	pH Buffer Calibration reagent
NIST Traceable:	Yes
Package Type:	Singlet
pH Value(s):	7.00 ± 0.02 pH at 25°C
Quantity:	20 pcs
Special Remarks:	Discard solution after use
Storage conditions:	10 - 25 °C (protect from light)