



Buffer Solution Kit, Colour-coded, pH 4.01, pH 7.00 and pH 10.01, 500 mL

Product #:

2947600

ZAR Price:

Contact Hach

Ships within 3 days

Save when you purchase pH buffer 3-packs!

pH Buffer Standard Solution for pH electrode calibration. Pack of 3 Colour-coded pH buffers. pH is ± 0.02 pH units. 500 mL bottles of pH 4.01, 7.00 & 10.01

33% Cost Reduction over single bottle purchase

Specifications

| | |
|---------------------|--|
| Color coded?: | Yes |
| Description: | pH Buffer Calibration reagent |
| NIST Traceable: | Yes |
| Package Type: | Bottle |
| pH Value(s): | 4.01, 7.00 & 10.01 ± 0.02 pH at 25°C |
| Quantity: | 3 pk of 500 mL |
| Storage conditions: | 10 - 25 °C (protect from light) |