



## Buffer Solution, pH 10.01, 500 mL

**Product #:** 1222149  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

### pH Buffer standard solutions for laboratory and field analyses

pH Buffer Standard Solution for pH electrode calibration. Colourless pH buffer solution. pH is  $10.01 \pm 0.02$  at  $25^{\circ}\text{C}$ . 500 mL bottle

#### Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

### Specifications

Color coded?:	No
Description:	pH Buffer Calibration reagent
NIST Traceable:	Yes
Package Type:	Bottle
pH Value(s):	$10.01 \pm 0.02$ pH at $25^{\circ}\text{C}$
Quantity:	500 mL
Special Remarks:	Observed shelf-life after opening, following recommended storage and use: 3 months
Storage conditions:	$10 - 25^{\circ}\text{C}$ (protect from light)