



## Buffer Solution Kit, Colour-coded, pH 4.01, pH 7.00 and pH 10.01, 20L

**Product #:** 2507300  
**ZAR Price:** Contact Hach  
No availability shown

### pH Buffer standard solutions for laboratory and field analyses

pH Buffer Standard Solution set for pH electrode calibration. Recommended 3-Pack of Colour-coded pH buffer solutions for easy identification. pH is  $\pm 0.02$  at 25°C. 20 L Bag-in-box packages of pH 4.01, 7.00 & 10.01, with a collapsible inner bladder in a white corrugated cardboard box. Ideal packaging size for laboratories doing high volume pH testing.

#### Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

#### Save when you purchase pH buffer 3-packs!

33% Cost Reduction over single bottle purchase

---

### Specifications

Color coded?:	Yes
Description:	pH Buffer Calibration reagent
NIST Traceable:	Yes
Package Type:	Bag-in-box (BiB)
pH Value(s):	4.01, 7.00 & 10.01 $\pm 0.02$ pH at 25°C
Quantity:	3 pk of 20 L
Storage conditions:	10 - 25 °C (protect from light)