



## Standard solution potassium hydroxide, 12.0 N, 500 mL

**Product #:**

**23049H**

ZAR Price:  
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For selenium determination by the Diaminobenzidine method

---

### Specifications

Concentration:	12.0 N
Platform:	Solution
Volume:	500 mL

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)