



Potassium hydroxide standard solution, 8.00 N, 500 mL

Product #: 28249
ZAR Price: Contact Hach
Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For the determination of Copper, 0.1-8.0 mg/L Cu, and Hardness, Titration and Buret Titration Method.

Specifications

Concentration:	8.00 N
Method:	EDTA Titration
Parameter:	varies
Platform:	Solution
Volume:	500 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)