



Sulphuric acid standard solution, 0.100 N, 1 L

Product #:

20253

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For Buret Acid-Base Titration, Chlorine Demand and Trihalomethane Formation Potential (THMFP).

Specifications

Concentration:	0.100 N
Method:	8288 (Acid)
Method Number 2:	8289 (Base)
Parameter:	varies
Platform:	Solution
Range:	10 - 100
Volume:	1 L

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)