



Sulfuric Acid Standard Solution, 1.000 N, 100 mL MDB

Product #:

127032

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Multi-purpose sulfuric acid standard solution for pH adjustment

Specifications

Application:	pH Adjustment
Description:	Buffer
Formula:	H ₂ SO ₄ 1 N
Package Type:	Bottle
Platform:	Solution
Quantity:	100 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)