



Sulfuric acid standard solution, 1.000 N, 1 L

Product #:

127053

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Multi-purpose sulfuric acid standard solution

Specifications

Concentration:	1.00 N
Description:	Sulfuric Acid Standard Solution
Method:	8288 (Acid)
Method Number 2:	8289 (Base)
Platform:	Solution
Range:	100 - 25000 meq/L
Volume:	1 L

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)