



Sodium hydroxide standard solution, 0.020 N (N/50), 1 L

Product #:	19353
ZAR Price:	Contact Hach
Available	

For acidity determination by buret titration. Hach methods 8219 (methyl orange) and 8010 (phenolphthalein). 1 liter bottle

Specifications

Concentration:	0.020 N
EPA compliant:	Yes
Method:	8219 (MO)
Method Number 2:	8010 (P)
Method Number 3:	8288 (Acid)
Method Number 4:	8289 (Base)
Platform:	Solution
Range:	0 - 5000 mg/L CaCO ₃ (MO)
Range 2:	0 - 10000 mg/L CaCO ₃ (P)
Range 3:	0 - 10 meq/L
Volume:	1 L

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)