



Sodium hydroxide standard solution, 0.100 N (N/10), 1 L

Product #: 19153
ZAR Price: Contact Hach
Ships within 3 days

For acidity determination by buret titration. 1 liter bottle

Specifications

Concentration:	0.100 N
Description:	Sodium Hydroxide Standard Solution, 0.100 N (N/10), 1 L
Method:	8288 (Acid)
Method Number 2:	8289 (Base)
Method Number 3:	8291
Platform:	Solution
Range:	10 - 100 meq/L
Volume:	1 L

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)