



Free Amino Nitrogen Reagent Set

Product #: 45950

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Simplify your fermentation efficiency by monitoring your Free Amino Nitrogen.

For use with your Hach DR6000 spectrophotometer and designed to simplify brewing analysis.

Testing for FAN is a part of standard brewing quality control analyses. This enables the protein content to be estimated and serves as an indicator of beer quality.

Specifications

Instrument:	DR6000
Parameter:	Free Amino Nitrogen
What's included?:	Di-Sodium-Hydrogen Phosphate, 500 g (170234)
	Potassium Phosphate, Monobasic, ACS, 454 g (17001H)
	Ninhydrin, 99%, 10 g (129523)
	D-Fructose, 99%, 250 g (135129)
	Potassium Iodate, ACS, 100 g (606932)
	Glycine, ACS, 100 g (211132)

What's included?

Di-Sodium-Hydrogen Phosphate, 500 g (170234)Potassium Phosphate, Monobasic, ACS, 454 g (17001H)Ninhydrin, 99%, 10 g (129523)D-Fructose, 99%, 250 g (135129)Potassium Iodate, ACS, 100 g (606932)Glycine, ACS, 100 g (211132)

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)