



## DR6000 VDK Hardware Kit

Product #: 4592510

**OBSOLETE ITEM**  
This item is no longer available.

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Avoid off-flavor in your beer and safeguard your product quality by measuring VDK. This package contains glassware you need to run the ASBC approved method.

---

### Specifications

Instrument:	DR6000
Parameter:	Vicinal Diketones
What's included?:	Package includes: 10 mL Volumetric Class A Flask; 50 mL Double Metric Cylinder; 50 mL Low Form Glass Beaker. (Reagent set ordered separately as PN 45910)

---

### What's included?

Package includes: 10 mL Volumetric Class A Flask; 50 mL Double Metric Cylinder; 50 mL Low Form Glass Beaker. (Reagent set ordered separately as PN 45910)