



DR6000 IBU Hardware Kit

Product #: 4592540

OBSOLETE ITEM
This item is no longer available.

Ensure your bitterness flavor profiles are precisely what you intend by measuring IBUs. This package contains everything you need to run the ASBC approved method.

Specifications

Instrument:	DR6000
Parameter:	Beer Bitterness
What's included?:	Package includes: Hydrochloric Acid; Iso-octane; 50 mL Disposable Centrifuge Tubes with Caps; 10 mm cuvette; 6 X 50 mL 115VAC Centrifuge; 120VAC Orbital Shaker; 500 mL Shaker Platform Flask Clamps; 1.0-10 mL Tensette Pipet; 1.0-10.0 mL Pipet Tips, 250/pk; 3.1 mL Bulb Draw Dropper, 20/pk

What's included?

Package includes: Hydrochloric Acid; Iso-octane; 50 mL Disposable Centrifuge Tubes with Caps; 10 mm cuvette; 6 X 50 mL 115VAC Centrifuge; 120VAC Orbital Shaker; 500 mL Shaker Platform Flask Clamps; 1.0-10 mL Tensette Pipet; 1.0-10.0 mL Pipet Tips, 250/pk; 3.1 mL Bulb Draw Dropper, 20/pk