



DR6000 Free Amino Nitrogen Hardware Kit

Product #: 4592550

OBSOLETE ITEM
This item is no longer available.

Simplify your fermentation efficiency by monitoring your Free Amino Nitrogen. This package contains lab apparatus you need to run the ASBC approved method. Reagents, Balance and pH meter and electrode sold separately.

Specifications

Instrument:	DR6000
Parameter:	Free Amino Nitrogen
What's included?:	Package includes: 10 mm cuvette; 100 mL Volumetric Class A Glass Flask; 118 mL Amber Glass Bottle, 6/pk; 22-400 foam liner caps; 1.0-10 mL Tensette Pipet; 1.0-10.0 mL Pipet Tips, 250/pk; Cimarec 7 x 7 hotplate; 125 mm X 65 mm Evaporating Glass Dish; Cell Holder Assembly; 10-20-25 mL Glass Sample Cell, 6/pk

What's included?

Package includes: 10 mm cuvette; 100 mL Volumetric Class A Glass Flask; 118 mL Amber Glass Bottle, 6/pk; 22-400 foam liner caps; 1.0-10 mL Tensette Pipet; 1.0-10.0 mL Pipet Tips, 250/pk; Cimarec 7 x 7 hotplate; 125 mm X 65 mm Evaporating Glass Dish; Cell Holder Assembly; 10-20-25 mL Glass Sample Cell, 6/pk

Required Accessories

- HQ411D Benchtop pH/mV Meter (Item HQ411D)
- Intellical PHC101 Laboratory Low Maintenance Gel Filled pH Electrode, 1 m Cable (Item PHC10101)
- ScienTech analytical balance - Zeta series - 120 g capacity (Item 2936801)