



DR6000 Total Polyphenols Hardware Kit

Product #: 4592570

OBSOLETE ITEM
This item is no longer available.

Hach confidence in your beer's color, taste and flavor profile by monitoring your Total Polyphenols. This lab package comes with hardware you need to run the ASBC approved method.

Specifications

Instrument:	DR6000
Parameter:	Total Polyphenols
What's included?:	Hardware Package includes:1000 mL, 100 mL, 25 mL Volumetric Class A Glass Flask; 0.1-1.0 mL & 1.0-10 mL Tensette Pipet; 0.1-1.0 mL & 1.0-10.0 mL Pipet Tips, 50/pk-250/pk
	Reagents available as Part Number 45970

What's included?

Hardware Package includes:1000 mL, 100 mL, 25 mL Volumetric Class A Glass Flask; 0.1-1.0 mL & 1.0-10 mL Tensette Pipet; 0.1-1.0 mL & 1.0-10.0 mL Pipet Tips, 50/pk-250/pkReagents available as Part Number 45970