



Reagent set acidity, digital titrator

Product #:

2272800

ZAR Price:

Contact Hach

Ships within 2 weeks



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Indicator Powder Pillows and Titration Cartridges for the determination of Acidity: Methyl Orange 10-160 mg/l CaCO_3 and Phenolphthalein 100-4000mg/l CaCO_3 .

For determination of acidity by Digital Titrator titration, with Methyl Orange (using bromphenol blue indicator) and Phenolphthalein endpoints. Hach Methods 8201 (Methyl Orange) and 8202 (Phenolphthalein). Ranges: 10-160 and 1000-4000 mg/L as CaCO_3 . Reagent Set includes: Bromphenol Blue Powder Pillows, Phenolphthalein Powder Pillows and two Sodium Hydroxide Cartridges (0.1600 N and 1.600N). Approximately 40 - 100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Method:	8201 (MO)
Method Number 2:	8202 (P)
Number of tests:	100
Platform:	Digital Titrator
Range:	10 - 160 mg/L
Range 2:	100 - 4000 mg/L CaCO_3

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)