



Sodium Hydroxide Digital Titrator Cartridge, 1.600 N

Product #:

1437901

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Ready-to-use titrant solution for Acid-Base Titration (1-10 meq/l acid) and determination of Acidity (100-4000mg/l CaCO_3).

For acidity determination by Digital Titration, 100 - 4000 mg/L CaCO_3

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

EPA compliant:	Yes
Method:	Fenolftaleína
Method Number 2:	8202 (P)
Method Number 3:	8200 (ácido)
Method Number 4:	8233 (Base)
Number of tests:	40 - 100
Parameter:	Acidity
Platform:	Digital Titrator
Range:	1 - 4000 mg/L CaCO_3
Range 2:	1 - 10 meq/L
Titrant:	1.600 N NaOH

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)