



Sodium Hydroxide Titration Cartridge 0.1600 N NaOH

Product #:	1437701
ZAR Price:	Contact Hach
Ships within 3 days	

For acidity determination by Digital Titration, 10 - 160 mg/L CaCO_3

Specifications

Method:	Methyl Orange
Method Number 2:	8202 (P)
Number of tests:	40 - 100
Parameter:	Acidity
Platform:	Digital Titrator
Range:	10 - 160 mg/L CaCO_3
Titrant:	0.1600 NaOH

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)