



Powder pillows indicator, bromphenol blue, pk/100

Product #: 1455099
ZAR Price: Contact Hach
Ships within 3 days

pH indicator for acidity and alkalinity titrations. Visual transition interval: pH 3.0 (yellow)-pH 4.6 (blue). Sample Size: 50 mL

Specifications

Method:	8219 (MO)
Method Number 2:	8010 (P)
Method Number 3:	8201 (MO)
Method Number 4:	8202 (P)
Number of tests:	100
Platform:	Powder Pillows
Range:	0 - 5000 CaCO ₃ (MO)
Range 2:	0 - 10000 CaCO ₃ (P)
Volume:	100 Powder Pillows