



Sulfuric acid cartridge for digital titrator, 0.1600 N

Product #: 1438801
ZAR Price: Contact Hach
Ships within 3 days

Ready-to-use titrant solution for the determination of alkalinity, 10 - 160 mg/L CaCO_3

Specifications

EPA compliant:	Yes
Method:	Phenolphthalein
Number of tests:	40 - 100
Parameter:	Alkalinity, Total
Platform:	Digital Titrator
Range:	10 - 160 mg/L CaCO_3
Titration:	0.1600 N H_2SO_4

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)