



Digital titrator cartridge, sulfuric acid, 1.600 N, 1 piece

Product #:

1438901

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Ready-to-use titrant solution for acid-base and alkalinity titrations.

Specifications

EPA compliant:	Yes
Method:	Phenolphthalein
Method Number 1 Chemistry:	8203
Number of tests:	40 - 100
Parameter:	Alcalinidad, base ácida
Platform:	Digital Titrator
Range:	100 - 4000 mg/L CaCO_3
Range 2:	1 - 10 meq/L Base
Titrant:	1.600 N H_2SO_4

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)