



Solution indicator, bromphenol blue, 100 mL MDB

Product #:

1455232

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

pH indicator for acidity and alkalinity titrations. Visual transition interval: pH 3.0 (yellow)-pH 4.6 (blue).

Specifications

Method:	Titration
Platform:	Solution
Volume:	100 mL MDB

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)