



Chlorophosphonazo Indicator Solution Pillows, pk/100

Product #:	2589599
ZAR Price:	Contact Hach
Available	

For the determination of Ultra Low Range Hardness, 8-1000 $\mu\text{g/L CaCO}_3$. Chlorophosphonazo indicator solution is packaged in unit-dose solution pillows.

For ultra low range hardness determination by the Chlorophosphonazo method. Hach Method 8374. Range: 4 - 1000 $\mu\text{g/L Ca \& Mg as CaCO}_3$. Chlorophosphonazo solution reagent packaged in solution pillows. Requires CDTA Reagent, Cat. No. 2589636. 100 Tests.

Specifications

Method:	Chlorophosphonazo
Method Number 1 Chemistry:	8374
Number of tests:	100
Parameter:	Calcium and Magnesium Hardness
Quantity:	100 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)