



Hardness Reagent Set, 8-1000 µg/L Ca & Mg as CaCO₃

Product #:

2603100

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For Ultra Low Range Hardness determination, Chlorophosphonazo Method.

Specifications

Instrument:	All except DR 900, DR/800 Series, PC II
Method:	8374
Number of tests:	100
Parameter:	Hardness
Platform:	Solution
Range:	8 - 1000 µg/L CaCO ₃
Special Features:	Reagent set, Pillows

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)