



Alkali solution for colorimetric hardness, 100 mL

Product #:

2241732

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For colorimetric Calcium and Magnesium determination, 0.05-4.00 mg/L Ca as CaCO₃, 0.05-4.00 mg/L Mg as CaCO₃.

Specifications

Method:	8030
Number of tests:	100
Parameter:	Hardness, Ca and Mg
Platform:	Solution
Quantity:	100 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)