



## Hardness Reagent Set, Calmagite Colorimetric Method

**Product #:**

**2319900**

ZAR Price:  
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For the determination of Hardness, 0.05-4.00mg/L Ca as  $\text{CaCO}_3$ /0.05-4.00mg/L Mg AS  $\text{CaCO}_3$ . 100 tests

---

### Specifications

Instrument:	All except PC II
Method:	8030
Number of tests:	100
Parameter:	Hardness, (Calcium and Magnesium)
Platform:	Solution
Range:	0.05 - 4.00 mg/L $\text{CaCO}_3$
Range 2:	0.05 - 4.00 mg/L Mg as $\text{MgCO}_3$
Special Features:	Reagent set

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)