



TitraVer EDTA Digital Titrator Cartridge, 0.800 M

Product #:	1439901
ZAR Price:	Contact Hach
Available	

Ready-to-use titrant solution for the determination of Calcium and Total Hardness, 100-4000 mg/L CaCO₃.

Specifications

Method:	EDTA (CalVer 2)
Method Number 2:	8213
Method Number 3:	8329
Number of tests:	40 - 100
Parameter:	Hardness, Total
Platform:	Digital Titrator
Range:	100 - 4000 mg/L CaCO ₃
Titration:	0.800 M EDTA

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)