



Hardness (Calcium) Reagent Set, LR, Digital Titrator

Product #:

2447200

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Reagent set for the determination of Calcium Hardness by EDTA Digital Titrator Titration Method (Hach method 8204). Range: 10-160 mg/L as CaCO_3 . Set contains CalVer 2 Calcium Indicator Powder Pillows, Potassium Hydroxide Solution (8 N) and 0.0800 M (EDTA) Digital Titrator Cartridge for approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Method:	EDTA Titration
Number of tests:	40-100
Parameter:	Hardness (Calcium)
Platform:	Digital Titrator
Range:	10 - 160 mg/L CaCO_3
Special Features:	Reagent set, LR, Digital Titrator

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)