



Total Hardness Reagent Set for Digital Titrator (LR)

Product #: 2448000
ZAR Price: Contact Hach
Ships within 3 days

For the determination of Total Hardness, 10-160 mg/L CaCO_3 , EDTA Method.

Specifications

EPA compliant:	Yes
Method:	EDTA Titration
Number of tests:	40-100
Parameter:	Hardness (Total)
Platform:	Digital Titrator
Range:	10 - 160 mg/L CaCO_3
Special Features:	Reagent set, LR, Digital Titrator

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)