



EDTA standard solution, 0.0200N, 4L

Product #:	20556
ZAR Price:	Contact Hach
Available	

For hardness determination by buret titration. Hach methods 8222 (calcium) and 8226 (total). 0.0200 N (0.0100 M) (N/50) 4 liter bottle

Specifications

Method:	EDTA Titration
Platform:	Solution
Quantity:	4
Range:	0 - 5,000 mg/L CaCO ₃

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)