



Hardness (Calcium & Total Sequential) HR Reagent Set, Digital Titrator

Product #:

2448700

ZAR Price:

Contact Hach

Ships within 3 days



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For sequential determination of Calcium and Total Hardness by EDTA Digital Titrator titration. Hach method 8329. Range: 100 - 4000 mg/L CaCO_3 . Set contains Buffer Solution (Hardness 1), CalVer 2 Calcium Indicator Powder Pillows, ManVer 2 Hardness Indicator Powder Pillows, Potassium Hydroxide Solution (8 N), Sulfuric Acid (5.25 N) and 0.800 M EDTA Digital Titrator Cartridge. Approximately 40 - 100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

| | |
|-------------------|---|
| Method: | 8329 |
| Number of tests: | 40 - 100 |
| Parameter: | Hardness (Calcium and Total Sequential) |
| Platform: | Digital Titrator |
| Range: | 100 - 4000 mg/L CaCO_3 |
| Special Features: | Reagent set, HR, Digital Titrator |

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)