



## Reagent set, hardness (total), buret, LR

<b>Product #:</b>	<b>2447600</b>
ZAR Price:	Contact Hach
Ships within 3 days	

Reagent set for the determination of Total Hardness by the ManVer 2 Buret Titration Method (Hach method 8226). EPA Accepted method. Range: 0-2,500 mg/L as CaCO<sub>3</sub>. Contains: Buffer Solution (Hardness 1), ManVer 2 Hardness Indicator Powder Pillows and 0.0200 N EDTA (TitraVer) Titrant for up to approximately 100 tests.

---

### Specifications

EPA compliant:	Yes
Method:	8226
Number of tests:	100
Parameter:	Hardness (Total)
Range:	0 - 2500 mg/L CaCO <sub>3</sub>
Special Features:	Reagent set, Buret, LR

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)