



HydraVer2 Hydrazine reagent solution, 4 - 600 µg/L, 100 mL

Product #:

179032

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For the determination of Hydrazine, p-Dimethylaminobenzaldehyde method. For 10 mL sample. 100 mL bottle.

Specifications

Description:	HydraVer 2 Hydrazine Reagent
Instrument:	All except DR900, PC II
Method:	p-Dimethylaminobenzaldehyde
Number of tests:	100
Parameter:	Hydrazine
Platform:	Solution
Quantity:	100 mL
Range:	4 - 600 µg/L N ₂ H ₄
Storage conditions:	10 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)