



## DEHA Reagent 1 Powder Pillows, 3-450µg/L DEHA

**Product #:** 2167969  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

Powder reagents in sealed packages for the determinaiton of oxygen scavengers by the iron reduction method. pk/100.

---

### Specifications

Description:	DEHA Reagent 1 Powder Pillows
Method:	Iron Reduction
Number of tests:	50
Parameter:	Oxygen Scavengers
Platform:	Powder Pillows
Quantity:	100/pk
Range:	3 - 450 µg/L DEHA
Range 2:	5 - 600 µg/L Carbohydrazide
Range 3:	9 - 1000 µg/L Hydroquinone
Range 5:	15 - 1000 µg/L Methylene ketoxime (MEKO)

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)