



Oxygen Scavenger Reagent Set, 3-450µg/L DEHA

Product #:

2446600

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Method available for a wide variety of Oxygen Scavengers: Carbohydrazide, DEHA, Hydroquinone, Iso-ascorbic (erythorbic) Acid, Methylethylketoxime. Ferrous Iron values are also determined. 100 tests

Specifications

Instrument:	All except DR/820, DR/850, and PC II
Method:	8140
Number of tests:	100
Parameter:	Oxygen Scavengers
Range:	5 - 600 Carbohydrazide
Range 2:	3 - 450 µg/L DEHA
Range 3:	9 - 1000 µg/L Hydroquinone
Range 4:	13 - 1500 µg/L Eritorbična kiselina / izoaskorbinska kiselina (ISA)
Range 5:	15 - 1000 µg/L Methylethylketoxime (MEKO)

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)