



Silica Reagent Set, 0.010-1.600 mg/L SiO₂

Product #:

2346800

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For the determination of Silica, Heteropoly Blue Method. For 25 mL sample. 100 tests

Specifications

Instrument:	DR/2000
Method:	Heteropoly Blue
Number of tests:	100
Parameter:	Silica
Range:	0.010 - 1.600 mg/L SiO ₂
Sample size:	25 mL
Special Features:	LR, 25 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)