



Reagent powder pillows, HR silica, acid, 10 mL, pk/1000

Product #: 2107428

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

Suggested replacements:
• 2107469

If necessary, we will change your selected shipping method to accomodate these items.

For determination of high range silica by the silicomolybdate method (Hach method 8185) Range:1 to 100 mg/L. Sample Size: 10 mL sample. Pack of 1000 Powder Pillows.

Specifications

Method:	Silicomolybdate
Number of tests:	1000
Platform:	Powder Pillows
Quantity:	1000 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)