



Silica Reagent Set, 1-100 mg/L SiO_2

Product #:

2429600

ZAR Price:
Available

Contact Hach



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Powder reagents in sealed packages for the determination of Silica, Silicomolybdate Method. For 10 mL sample: pk/100.

Specifications

Instrument:	All except DR/3000, DR/2000, DR/850, DR/820, PC II
Method:	Silicomolybdate
Number of tests:	100
Parameter:	Silica
Platform:	Powder Pillows
Range:	1 - 100 mg/L SiO_2
Sample size:	10 mL
Special Features:	HR, 10 mL, 100 tests

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)