



Silica Reagent Set, 0.010-1.600 mg/L SiO₂

Product #:

2459300

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For the determination of Silica, Heteropoly Blue Method. For 10 mL sample.

Specifications

Instrument:	All except DR/3000, DR/2000, PC II
Method:	Heteropoly Blue
Number of tests:	100
Parameter:	Silica
Range:	0.010 - 1.600 mg/L SiO ₂
Sample size:	10 mL
Special Features:	LR, 10 mL, 100 tests

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)