



Acid reagent powder pillows, 25 mL, pk/1000

Product #:

104228

ZAR Price:

Contact Hach

No availability shown

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For determination of high range silica by the silicomolybdate method (Hach method 8185) Range: 1-100.0 mg/L SiO₂. Also used in the iodometric method (Hach method 10100) for the determination of high range hypchlorite (50-150 g/L as Cl₂). Pack of 1000 powder pillows.

Specifications

Method:	Silicomolybdate
Platform:	Powder Pillows

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)