



Reagent set silica, low range, 25 mL

Product #:

2256900

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For low range silica determination by the Heteropoly Blue method. 25 mL sample Approximately 100 tests

Specifications

Method:	Heteropoly Blue
Number of tests:	100
Range:	0.010 - 1.600 mg/L SiO_2 and 3 - 1,000 $\mu\text{g/L}$ SiO_2
Sample size:	25 mL
Special Features:	LR, 25 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)