



Reagent powder pillows, silica, citric acid, 10 mL, pk/1000

Product #:	2106228
ZAR Price:	Contact Hach
Ships within 3 days	

For determination of high range silica by the silicomolybdate method (Hach method 8185), and low range silica by the heteropoly blue method (Hach method 8186). Pack of 1000 Powder Pillows.

Specifications

Parameter:	Silica
Platform:	Powder Pillows
Quantity:	pkg 1000

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)