



## Molybdate 3 reagent solution, 100 mL

**Product #:**

**199532**

**ZAR Price:**

Contact Hach

**Ships within 3 days**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For the determination of low range (0 - 1.60 mg/L SiO<sub>2</sub>) and ultra low range (0 - 1000 ug/L SiO<sub>2</sub>) silica, Heteropoly Blue Method.

---

### Specifications

Method:	Heteropoly Blue
Method Number 2 Chemistry:	8282
Number of tests:	100
Parameter:	Silica
Platform:	Solution
Quantity:	100 mL

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)