



Citric acid reagent solution for rapid liquid Silica, 500 mL

Product #:

2254249

ZAR Price:
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For the determination of ultra low range Silica, 3 - 1000 µg/L SiO₂, Heteropoly Blue Method.

Specifications

Method:	8282
Number of tests:	250
Parameter:	Silica
Platform:	Solution
Quantity:	500 mL

Required Accessories

- Reagent powder, amino acid F, rapid liquid silica, 55 g (Item 2651155)
- Reagent solution, molybdate 3, 500 mL (Item 199549)
- Amino acid F dilution solvent, for rapid liquid silica, 475 mL (Item 2353011)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)