



Be Right™



## Silica Reagent Set, Bulk Solution, 3-1000 µg/L SiO<sub>2</sub>

Product #:

2553500

ZAR Price:

Contact Hach

Available

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

For the determination of ultra low range Silica, Heteropoly Blue Method. Approximately 100 tests.

---

### Specifications

Instrument:	All except DR 900, DR/800 Series, PC II
Method:	8282
Number of tests:	100
Parameter:	Silica
Platform:	Solution
Range:	3 - 1000 µg/L SiO <sub>2</sub>

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)